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Safety valve NO PVC

CE 1014

Catalogue No.:	V696433
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designed user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	ABS, rubber
Compatibility:	With medical devices equipped with standard Female Luer Lock connector in conformity with standard ČSN EN ISO 80369-7. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

The safety valve is designed for intermittent injection of drugs through a rubber latex-free stopper that allows repeated application and protects against bacterial contamination. Safety valves can be used for both gravity and pressure infusion systems up to a pressure of 2 bar.

Instruction for use

- Under aseptic conditions remove the safety valve from the packaging.
- Before application wipe or spray the surface of the rubber stopper with an alcoholic disinfectant.
- Wait for no less than 30 seconds until the disinfectant solution dries.
- Attach the valve with a clockwise twisting motion so as to ensure a safe connection.
- Detach the valve with a counter clockwise twisting motion.
- After use, products should be disposed of in line with prescribed procedures for disposal of medical (infectious) waste in accordance with local regulations.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain PVC
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Pressure
	Medical Device		Date of Manufacture				

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

 GAMA GROUP a.s., Betrieb Jimramov – 592 42 Jimramov
Sitz der Gesellschaft GAMA GROUP a.s., Mánesova 11/3b, České Budějovice 7, 370 01 České Budějovice, Tschechische Republik, Id.Nr.: 45797463, www.gama.cz

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