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Urologic catheter for women Catheter Nelaton Catheter Tiemann

Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: population of pediatric or adult patients
Designated user: medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used: PVC
Compatibility: With universal gradual (step) connector. With irrigation (lavage) syringe Janette. The catheters are not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the catheter is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date.
 Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

Catheters are intended to be introduced via ureter into the bladder cavity and to void and/or rinse the bladder.

Instruction for use

1. Check the type labelling of the catheter, usable life on the package and the package for damages.
2. Take the catheter out of the packaging under aseptic condition.
3. Use the catheter for the intended purpose.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Mountain	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen free	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Do not Resterilize	Do not Resterilize	Single use only	For single use only	Do not use if unit package is damaged	Do not use if unit package is damaged	MD	Medical Device
Date of Manufacture	Date of Manufacture						

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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