CS Combi zátky LL

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EN LL combo stoppers

Catalogue No.: Z676004, Z676005, Z676006
Medical device variants: LL combo stopper – red NO PVC, LL combo stopper – blue NO PVC, LL combo stopper – white NO PVC
Product type: medical device
Design: sterile, ethylene oxide sterilization
Target patient: for all patients, without restrictions regarding gender, age and can be used for adults, children and newborns
Designated user: healthcare professional
Materials used: PP
Compatibility: With medical devices equipped with standard Luer Lock connectors in conformity with standard ČSN EN ISO 80369-7. LL combo stopper is not intended for drug delivery.
Contraindication: Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects: No adverse effects are known if the combo stopper is used in accordance with the instructions specified in the Warning and Use sections.
Warning: Do not use, if the sterile packing is opened or damaged. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage: The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

Single use medical devices designed for safety closure where a Luer closure is required. The combination of MLL and FLL connection cones allows safe and versatile use when closing syringes, intravenous catheters, three-way cocks, ramps and other infusion systems. Combi stoppers are designed for gravity systems.

Instruction for use

- Under aseptic conditions remove the stopper from the packaging.
- Attach the stopper with a clockwise twisting motion so as to ensure a safe connection.
- Detach the stopper with a counter clockwise twisting motion.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	Date	Date of Manufacture
Resterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	MD	Medical Device

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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