CS

GAMAPLUS 3,0 x 150 LCH / GAMAPLUS 4,5 x 20 LCH

CE 1014

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EN GAMAPLUS 3,0 x 150 LCH / GAMAPLUS 4,5 x 20 LCH

CE 1014

| | |
|-------------------|--|
| Catalogue No.: | V606372-ND / V606372-01ND |
| Product type: | medical device |
| Design: | sterile, ethylene oxide sterilization |
| Target patient: | population of pediatric or adult patients |
| Designated user: | medical personnel with appropriate education, knowledge and practical experience in medical devices. |
| Materials used: | PVC, PE |
| Compatibility: | With universal gradual (step) connectors. With Female Luer Lock connector which is in accordance with the standard ČSN EN ISO 80369-7. GAMAPLUS is designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see www.sukl.cz or www.olecich.cz |
| Contraindication: | Do not use on any patient with a known allergic reaction to any of the product components. |
| Adverse effects: | No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections. |
| Warning: | Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient. |
| Storage: | The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. |

Intended use

Gamplus is designed to connect two different transfer sets which are ended with the negative cone and the gradual connecting cone.

Instruction for use

- Under aseptic conditions, remove the item from the packing.
- Remove the positive cone cover.
- Connect the positive cone to the negative cone of the transfer set.
- Connect the white terminal to the gradual connecting cone of the transfer set.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

| REF | Catalogue Number | Keep away from Rain | LOT | Batch Number | Manufacturer |
|-----|----------------------|---------------------|-----|---------------------------------------|---|
| | Use by / Expiry date | | | Read instructions for use | STERILE EO Sterilised by Ethylene Oxide |
| | Pyrogen free | | | Latex free product | Does not contain DEHP phthalate |
| | Do not Resterilize | | | Do not use if unit package is damaged | Pressure |
| | Medical Device | | | Date of Manufacture | Volume of tubing |

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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