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## Infusion set IS-121 / Infusion set IS-123

<b>Catalogue No.:</b> V606121 V606123	<b>Medical device variants:</b> Infusion set IS-121 Infusion set IS-123
<b>Product type:</b>	medical device
<b>Design:</b>	sterile, ethylene oxide sterilization
<b>Target patient:</b>	population of pediatric or adult patients
<b>Designated user:</b>	medical personnel with appropriate education, knowledge and practical experience in medical devices.
<b>Materials used:</b>	PVC, PE, PS, PP, ABS, SILICONE
<b>Compatibility:</b>	With medical devices equipped with standard Luer connectors in conformity with standard ČSN EN ISO 80369-7. With plastic infusion bags, plastic or glass bottles. With infusion pharmaceuticals without any special requirements for filtration. (In)compatibility of the drugs with the materials used check in the summary about products (SPC) see <a href="http://www.sukl.cz">www.sukl.cz</a> or <a href="http://www.olecich.cz">www.olecich.cz</a> With infusion pharmaceuticals incompatible with PVC (cytostatics) and with infusion pharmaceuticals sensitive to light.
<b>Incompatibility:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components.
<b>Adverse effects:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
<b>Warning:</b>	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Infusion bottle of contact with the infusion set until immediately before administration. During the operation, check the patient's condition and functionality sets. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
<b>Storage:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

<b>Intended use</b>	The sets are envisaged to administer of infusion from bottles or bags using infusion pumps.	Set:	IS-121 IS-123	Infusion pump:	IVAC Star-flow 581 IVAC Star-flow 591
---------------------	---	------	------------------	----------------	--

**When using the sets see to it that the attached pumps are of the stipulated type.**

### Instruction for use

- Under aseptic conditions, remove the set from the packing.
- Close the plug in the aeration hole on the body of the piercing point.
- Close the control clamp on the hose downstream of the drip chamber.
- Detach the cap from the piercing point and remove the cone cover.
- Insert the piercing point through the stopper to a bottle or a bag. Hang the bottle or the bag on a rack.
- Press the drip chamber so that it fills up to one half, open the control clamp, open the plug in the aeration hole and drive air out of the hose.  
Note: In case the set is to transport saline from the plastic bag, keep the aeration hole closed.
- Close the control clamp.
- Then follow the instructions on use of the pump to be installed jointly with the set.

**NOTE:** With the pump operational, the control clamp must be fully open. With the dosing process finished, the control clamp is to be closed before the set may be removed from the device.

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

	Catalogue Number		Keep away from Rain		Batch Number		Manufacturer
	Use by / Expiry date		Temperature Limitation		Read instructions for use		Sterilised by Ethylene Oxide
	Pyrogen free		Keep away from sunlight		Latex free product		Does not contain DEHP phthalate
	Do not Resterilize		For single use only		Do not use if unit package is damaged		Pressure
	The number of drops per milliliter		Filter of liquid with pore size		Medical Device		Date of Manufacture
	Volume of tubing		Contains phthalate DEHP				

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

### INFORMATION FOR THE RISK GROUP OF PATIENTS (This notice applies to product V606121)

Should treatment using this medicinal product involve treatment of children, pregnant or breast-feeding women, please become familiar with the following possible residual risks which may arise if using the product.

#### Residual risks with the application of this medicinal product

The product contains phthalates. If phthalates are inhaled, consumed or penetrate the skin, they may cause cancer or increase its incidence. They may also cause congenital genetic damage or increase its incidence and may cause or increase the incidence of non-congenital negative effects on offspring or exacerbate the male or female reproductive functions or abilities.

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