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## Venous pressure meter

<b>Catalogue No.:</b> V646992-ND	<b>Medical device variants:</b> Venous pressure meter MŽT-01	<b>Catalogue No.:</b> V646992-01ND	<b>Medical device variants:</b> Venous pressure meter MŽT-01 MLL	<b>Catalogue No.:</b> V646993-ND	<b>Medical device variants:</b> Venous pressure meter MŽT-02
<b>Functionality:</b>	the manufacturer declares the functionality of the device when used in accordance with the instructions for use				
<b>Duration of use:</b>	„Short term“ - continuous use for the time maximum is 24 hours				
<b>Design:</b>	sterile, ethylene oxide sterilization				
<b>Target patient:</b>	population of pediatric or adult patients				
<b>Designated user:</b>	a healthcare worker with adequate education, knowledge and practical experience focused on this MD. The product is intended exclusively for professional use.				
<b>Materials used:</b>	PE, PVC, PC – all types, ABS, elastomer – only MŽT-02				
<b>Compatibility:</b>	With infusion sets in conformity with standard ČSN EN ISO 8536-4 or ČSN EN ISO 8536-8. With intravenous catheters equipped with a connector with Luer thread in conformity with standard ČSN EN ISO 80369-7. The venous pressure meter is not intended for drug delivery.				
<b>Contraindication:</b>	Do not use on any patient with a known allergic reaction to any of the product components. Measurements under various conditions (change of position adjustment zero).				
<b>Adverse effects:</b>	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.				
<b>Warning:</b>	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Reuse, re-sterilization, or non-observance of storage conditions creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.				
<b>Storage and transport:</b>	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%. Organic solvents and chemicals must not be stored with the products. If the package is exposed to conditions outside the specified conditions, contact the manufacturer who will decide on the next course of action. The products are transported in clean, dry and covered means of transport at a temperature from -15°C to +60°C and a relative humidity of max. 75% so that they are not damaged. No chemicals may be transported with them. The maximum transport time of MD at temperatures of +50°C to +60°C must not exceed 4 days!				

### Intended use

The venous pressure meter is a single-use medical device used to orientatively determine the patient's central venous pressure (CVP). The value of the patient's central venous pressure indirectly informs about the level of hydration and the size of venous return and is therefore only indicative. This is not a treatment procedure, but a way of monitoring the trend of CVP. The scale on the manometric branch is only indicative and is not used as a gauge it only serves to determine the approximate value of CVP. A central venous catheter is used for accurate measurement of the CVP value.

The venous pressure meter may be applied only when accessorized with an auxiliary device (not delivered as a part of the product) capable of providing:

- Vertical position of the related pressure-gage branch with a scale printed on it;
- Adjustment of scale zero point (0) so that the point corresponds with the height of jugular depression (fossa jugularis) in a lying patient.


### Instruction for use

- Take the venous pressure meter out of the packing under aseptic conditions and place it onto the auxiliary device (see Figure 1).
- Remove the cover from the three-way valve and attach the branch tube to the infusion set. Take off the hose cover; the hose terminates in:
  - a positive cone (type: MŽT-01), or
  - a "Y" piece with a rubber stopper and a luer lock cone (type: MŽT-02).
- Turn the three-way valve to Position 2 (see Figure 3) and keep it there long enough for the pressure gage branch of the venous pressure meter to be filled up with the solution to the level of approx. 10 cm above the zero point of the scale (see Figure 1).
- Turn the three-way valve to Position 3 (see Figure 4) and read the solution level on the scale of the pressure gage branch (in cm). The value of the level of the solution corresponds to the approximate value of the venous pressure in the patient's central veins..

### PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

### Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use		Sterilised by Ethylene Oxide and single sterile barrier system
Pyrogen free		Keep away from sunlight	Latex free product		Does not contain DEHP phthalate
Do not Resterilize		For single use only	Do not use if unit package is damaged		Medical Device
Date and Country of Manufacture		UDI	Unique Device Identifier		

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