## Rampa s bezjehlovým vstupem

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## Ramp with needleless port

**Catalogue No.:** V696436  
**Product type:** medical device  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** population of pediatric or adult patients  
**Designated user:** medical personnel with appropriate education, knowledge and practical experience in medical devices.  
**Materials used:** PSU, PC, PE  
**Compatibility:** With medical devices equipped with standard Luer connectors under ČSN EN ISO 80369-7. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz)  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the ramp is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Prior to use, check all the taps for operability. Use standard Luer connection only, which complies with harmonised standards, so as not to damage the needleless ports. When using the ramps it is necessary to adhere to the principles of antiseptics, sterility and the specified replacement period must be respected.  
 Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date.  
 Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
 If the ramp is not exchanged after 72 hours of use, cracking may occur resulting in air embolism, leaks of the infusion solution and/or its contamination..  
**Storage:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.  
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

**Intended use**  
 The needleless port ramps (hereinafter referred to as the Ramp) are designed for simultaneous infusion therapies to administer infusions or drugs. The ramp is designed for connection to the central venous catheter for pressure and gravity applications at the pressure of up to 1.5 bar. It allows configuration, management, control or suspension of the fluid's desired flow direction. The needleless connection cones are compatible with Luer cones. The ramp must be always replaced no later than after 7 days or 100 connections (one needleless port = 100 connections).

**Instruction for use**  
 1. Under aseptic conditions remove the ramp from the packaging and remove the guard from the connection cone.  
 2. Before application, wipe or spray the surface of the silicone diaphragm of needleless ports with an alcoholic disinfectant.  
 3. Wait for no less than 30 seconds until the disinfectant solution dries out.  
 4. Connect the ramp using the connection - extension tubing (not included in the kit) to a central venous catheter to ensure a safe connection.  
 5. Attach the needleless ports to infusion or injection lines with Luer ending so as to ensure a safe connection.  
 6. Turn the taps to set the desired direction of flow.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**  
 Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

| Used symbols on packaging |                      |  |                         |  |                                       |  |                              |
|---------------------------|----------------------|--|-------------------------|--|---------------------------------------|--|------------------------------|
|                           | Catalogue Number     |  | Keep away from Rain     |  | Batch Number                          |  | Manufacturer                 |
|                           | Use by / Expiry date |  | Temperature Limitation  |  | Read instructions for use             |  | Sterilised by Ethylene Oxide |
|                           | Pyrogen free         |  | Keep away from sunlight |  | Latex free product                    |  | Does not contain PVC         |
|                           | Do not Resterilize   |  | For single use only     |  | Do not use if unit package is damaged |  | Pressure                     |
|                           | Medical Device       |  | Date of Manufacture     |  |                                       |  |                              |

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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