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Flow reversal device

Catalogue No.:	V606136-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PE, PVC, PP
Compatibility:	With dialysis sets in conformity with standard ČSN EN ISO 8637-2. With end pieces of vascular access having threads in conformity with standard ČSN EN ISO 80369-7. The Flow reversal device is not intended for drug delivery.
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Description of the method

According to specialised guidelines (K/DOQI¹, EBPG²), regular assessment of blood flow (QVA) is an optimal method of monitoring the quality of vascular access for the purpose of haemodialysis. Majority of QVA-assessment methods requires short-time inverted connection of needles to an extracorporeal circuit (an arterial needle connected to a venous blood set, and a venous needle connected to an arterial set). The Flow Reversal Device allows inverting the connection of needles to an extracorporeal circuit by means of a simple change of closed clamps while no disconnection is required (refer to Fig. 1). Thus, it considerably facilitates the measurement process, and increases the hygiene of work. The set is installed between the needles and respective outlets from the blood set (red connectors of the Flow Reversal Device are connected to the arterial needle and the arterial set; blue connectors to the venous needle and the venous set), and can be used in combination with any type of a blood set. In case of normal connection, the red arterial clamp (A) and the blue venous clamp (V) are open during dialysis; in inverted connection, A, V clamps are closed, and inverting white RA, RV clamps are open.

The Flow Reversal Device is a medical device of IIa class.

¹) K/DOQI Guidelines on Vascular Access, Am J Kidney Diseases, 48, Suppl 1, 2006

²) European Best Practice Guidelines on Vascular Access, Nephrol Dial Transpl, 22, 2007, Suppl 2

Intended use

Single use medical device intended for reversion the connection of needles in vascular access to an extracorporeal circuit upon vascular access blood flow measurement during hemodialysis.

Instruction for use

1. The preparation of an extracorporeal circuit with a Flow Reversal Device

Under aseptic conditions, remove the set from the packing. Prior to flushing the arterial set with physiologic or substitution solution, connect the red connector of the Flow Reversal Device to the respective connector of the arterial blood set. Open the solution inlet to the arterial set, and fill the entire Flow Reversal Device with the solution by gradually opening individual clamps. Then flush the entire extracorporeal circuit with a dialyser normally.

2. Commencement of dialysis

Prior to the commencement of dialysis, check the closure of all clamps on the Flow Reversal Device and complete filling of all its parts with the solution. Connect the red connector of the Flow Reversal Device to the arterial needle, open the red A clamp of the Flow Reversal Device and the clamp on the inlet part of the standard arterial set and on the outlet part of the venous set, and start the blood pump. After filling up the extracorporeal circuit with blood up to the optical detector of the dialysis apparatus, stop the blood pump, connect the venous set to the respective blue connector of the Flow Reversal Device and the other blue connector of the Flow Reversal Device to the venous needle. Open the blue V clamp of the Flow Reversal Device and the clamp on the venous set, and re-start the blood pump.

Warning:

After the commencement of dialysis, duly fix the Flow Reversal Device with adhesive plaster at the connection of needles to the forearm, or to an arm, so that the withdrawal of the needles due to the weight of Flow Reversal Device is prevented!!

3. Switching normal connection to inverted connection and vice versa during dialysis

When switching from normal to inverted connection upon vascular access blood flow measurement during dialysis, first, open both white RA, RV clamps, and immediately afterwards close both coloured clamps (the red A clamp and the blue V clamp). During this manoeuvre, neither the blood pump has to be stopped (some systems of vascular access blood flow measurement do not even allow this), nor the set blood flow in the extracorporeal circuit has to be reduced. If the switching is quick enough, the alarms of pressure sensors of the extracorporeal circuit are not usually activated. When switching from inverted connection back to normal connection, the process is similar – first, open both coloured clamps, and then close both the white inverting clamps of the Flow Reversal Device.

Warning:

Theoretically, the Flow Reversal Device allows 16 different combinations of open and closed clamps, out of which only two are utilised (A, V open and RA, RV closed – normal connection; A, V closed, RA, RV open – inverted connection). Incorrect combination of open and closed clamps (improper closure of a clamp) always causes certain change in pressure in the extracorporeal circuit. Such a change, however, may not always be sufficient to activate a pressure sensor alarm. Nevertheless, in such case increased recirculation of blood in the extracorporeal circuit may occur resulting in the decrease in dialysis efficiency and the risk of subsequent blood coagulation. Thus, when switching clamps, always pay maximum attention to pressure changes in order to detect such fault. If the dialysis apparatus is equipped with a pressure sensor between the blood pump and the dialyser, always make use of it. It increases the possibility of detection of wrong or improper opening/closing of a clamps of the Flow Reversal Device!!

Note 1:

When using ultrasonic device Transonic HD01/ 02 (generally considered as reference measurement of vascular access blood flow), the device ultrasonic sensors need to be fitted on the tubes of the Flow Reversal Device behind crossing, i.e. closer to the connectors connecting the Flow Reversal Device to the arterial and venous blood sets, not on the part before crossing, i.e. closer to the connection of the set to needles !!

Note 2:

Provided that you make use of a device only capable of measuring, in inverted connection, forced recirculation, but not the blood flow proper, the vascular access blood flow value may be determined based on the value of the forced recirculation measured and the extracorporeal circuit blood flow by means of the graph in Fig. 2.

4. Dialysis termination

Dialysis shall always be finished with open A, V clamps and closed RA, RV clamps. If blood from the arterial set is to be returned by air, always disconnect the connector between the Flow Reversal Device and the needle, not between the Flow Reversal Device and the blood set in order to minimize the blood loss. After the return of blood, stop the blood pump, close the blue V clamp of the Flow Reversal Device, and disconnect the Flow Reversal Device from the venous needle.

Note 3:

Procedure for safe disposal of medical devices after use:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.


Fig. 1 Flow Reversal Device


Fig. 2

Used symbols on packaging

REF	Catalogue Number	Keep away from Rain	LOT	Batch Number	Manufacturer
Use by / Expiry date		Temperature Limitation	Read instructions for use	Sterilised by Ethylene Oxide	
Pyrogen free		Keep away from sunlight	Latex free product	Does not contain DEHP phthalate	
Do not Resterilize		For single use only	Do not use if unit package is damaged	Medical Device	
Date of Manufacture					

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¹⁾ K/DOQI Guidelines on Vascular Access, Am J Kidney Diseases, 48, Suppl 1, 2006
²⁾ European Best Practice Guidelines on Vascular Access, Nephrol Dial Transpl, 22, 2007, Suppl 2

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