## Spojovací hadička

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**EN Connecting tubing**

**Catalogue No.:** V606402-ND  
 V606404-ND  
**Medical device variants:** Connecting Tubing 2-way  
 Connecting Tubing 4-way  
**Product type:** medical device  
**Design:** sterile, ethylene oxide sterilization  
**Target patient:** population of pediatric or adult patients  
**Designated user:** medical personnel with appropriate education, knowledge and practical experience in medical devices.  
**Materials used:** PVC, PE  
**Compatibility:** With medical devices equipped with standard Luer connectors in conformity with standards ČSN EN ISO 80369-7. Connecting Tubing are designed to transfer liquid or i.v. administration or transfer of blood. (In)compatibility of the drugs with the materials used, check in the summary about products (SPC) see [www.sukl.cz](http://www.sukl.cz) or [www.olecich.cz](http://www.olecich.cz).  
**Contraindication:** Do not use on any patient with a known allergic reaction to any of the product components.  
**Adverse effects:** No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.  
**Warning:** Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date.  
 Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.  
**Storage:** The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation.  
 Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

**Intended use**

Connecting tubing is used to connect more infusion sets with other components, such as intravenous cannula. Inputs of connecting tubes made up by Luer Lock connecting cones can be used to administer one or more infusion or drug solutions. The connecting tubes are designed for both pressure and gravity systems up to a pressure of 2 bar and are compatible with Luer-Lock cones, which meet the requirements of ČSN EN ISO 80369-7.

**Instruction for use**

- Under aseptic conditions remove the tube from the packaging.
- Remove the caps from the connecting beams.
- Attach the tube's negative cone to the infusion set (or syringe) and the positive cone to the intravenous cannula.

**PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:**

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

**NOTE TO USERS:** Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

**Used symbols on packaging**

REF	Catalogue Number	Umbrella	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Resterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	P	Pressure
MD	Medical Device	Date of manufacture	Date of Manufacture	VOL	Volume of tubing		

The instruction for use is also available in electronic form at [www.gama.cz](http://www.gama.cz). Free copies of other copies are available on request via e-mail: [sales@gama.cz](mailto:sales@gama.cz).

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