Urodynamická souprava malá USM

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Small urodynamic set USM

Catalogue No.:	V606481-ND
Product type:	medical device
Design:	sterile, ethylene oxide sterilization
Target patient:	population of pediatric or adult patients
Designated user:	medical personnel with appropriate education, knowledge and practical experience in medical devices.
Materials used:	PVC, PE, PP, PS, ABS
Compatibility:	With urodynamic appliances UROMIC series. With cystometry urological catheter. With plastic infusion bags, plastic or glass bottles. Small urodynamic set USM is not intended for drug delivery.
Contraindication:	Do not use on any patient with a known allergic reaction to any of the product components.
Adverse effects:	No adverse effects are known if the set is used in accordance with the instructions specified in the Warning and Use sections.
Warning:	Do not use if the sterile packing is opened or damaged, or when the protective covers dropped off. Do not use after the expiry date. Repeated use or repeated sterilization creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability, may lead to injury, illness or death of the patient.
Storage:	The products should be stored in dry, ventilated, dust-free and dark rooms without direct sunlight to avoid their damage and deformation. Products must be stored in the packaging as delivered by the producer. Products must be stored at temperatures from +5°C to +25°C and the relative air humidity is max. 75%.

Intended use

Single use medical device intended as accessory to UROMIC apparatus during urodynamic measurements.

Use the sets only in connection with the prescribed type of pump.

Instruction for use

- Under aseptic conditions, remove the set from the packing. Close the pinchcock on the tube.
- Remove the cover from the puncturing style and the cover from the cone.
- Insert the puncturing style vertically into the stopper of the bottle or sac. Hand the bottle or sac on the stand.
- Fill up the drip chamber halfway.
- Open the stopper of the aeration opening, open also the pinchcock - and remove any air from the hose.

Note: If the set is used to transfer solution from the plastic pack, leave the aerating inlet closed.

- Close the pinchcock on the tube.
- Place the flexible part of the set into the pump, open the pinchcock.
- Connect the attachment cone to the opening of the catheter.
- Follow the instructions for using the pump for which the kit is designated.

NOTE: When the pump is running, the pinchcock must be fully open. After the dosing is finished, it is necessary to close the regulating pinchcock before the set is removed from the apparatus.

PROCEDURE FOR SAFE DISPOSAL OF MEDICAL DEVICES AFTER USE:

Put the medical device in a blue PE bag intended for hazardous - infectious waste. After filling, seal the labeled bag (waste code, place of origin, date and signature of the responsible person) and transport it safely to the incinerator = decontamination device. Removal of the dangerous property of infectivity means proper decontamination. Waste sorting is carried out according to the operating rules of the medical facility.

NOTE TO USERS: Any serious adverse event occurring in relation to the device in question should be reported to the manufacturer and the competent local authority.

Used symbols on packaging

REF	Catalogue Number	Rain	Keep away from Rain	LOT	Batch Number	Factory	Manufacturer
Hourglass	Use by / Expiry date	Thermometer	Temperature Limitation	Book	Read instructions for use	STERILE EO	Sterilised by Ethylene Oxide
Pyrogen	Pyrogen free	Sun	Keep away from sunlight	Latex	Latex free product	DEHP	Does not contain DEHP phthalate
Resterilize	Do not Resterilize	Single use	For single use only	Damaged	Do not use if unit package is damaged	Pressure	Pressure
MD	Medical Device	Date of Manufacture	Date of Manufacture				

The instruction for use is also available in electronic form at www.gama.cz. Free copies of other copies are available on request via e-mail: sales@gama.cz.

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